


Product Name: Recombinant Human ARID5B Protein
Catalog #: PEH4915



Summary

Name	ARI5B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14865
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTPQGRNSDH GEDEVIADVSE KVIVKLEDLV KVVHSDFSKW RCGFHAGPVK TEALGRNGQK EALLKYRQST LNSGLNFKDV LKEKADLGED

image



Background

Alternative Names AT-rich interactive domain-containing protein 5B (ARID domain-containing protein 5B) (MRF1-like protein) (Modulator recognition factor 2) (MRF-2)

Note

For Research Use Only , Not for Diagnostic Use