


Product Name: Recombinant Human ATXN7L3 Protein
Catalog #: PEH4917



Summary

Name	AT7L3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14CW9
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLEKCLGMGR NSSRIANRRI ANSNMNMKSE SDQEDNDDIN DNDWSYGSEK KAKKRKSDKN PNSPRRSKSL KHKNGELSNS DPFKYNNSTG ISYETLGPEE LRSLTTQCG VISEHTKKMC TRSLRCPQHT DEQRRTVRIY FLGPSAVLPE VESSLDNDSF DMTDSQALIS RLQWDGSSDL SPDSGSSKT SENQGWGLGT NSSESRTKK KSHLSLVGT ASGLGSNKKK KPKPPAPPTP

image



Background

Alternative Names	Ataxin-7-like protein 3 (SAGA-associated factor 11 homolog)
--------------------------	---

Note

Product Name: Recombinant Human ATXN7L3 Protein
Catalog #: PEH4917



For Research Use Only , Not for Diagnostic Use