



Product Name: Recombinant Human TAF1C Protein
Catalog #: PEH4921



Summary

Name	TAF1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15572
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEPGRQTRRP KRRTQLSSSF SLSGHVDPSE DTSSPHSPEW PPADALPLPP TTPPSQELTP DACAQGVPS E QRQMLRDYMA KLPPQRDTPG CATTPPHSQA SSVRATRSQQ HTPVLSSSQP

image



Background

Alternative Names	TATA box-binding protein-associated factor RNA polymerase I subunit C (RNA polymerase I-specific TBP-associated factor 110 kDa) (TAFI110) (TATA box-binding protein-associated factor 1C) (TBP-associated factor 1C) (Transcription initiation factor SL1/TIF-IB subunit C)
--------------------------	---

Note

Product Name: Recombinant Human TAF1C Protein
Catalog #: PEH4921



For Research Use Only , Not for Diagnostic Use