


Product Name: Recombinant Human MED1 Protein
Catalog #: PEH4922



Summary

Name	MED1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15648
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FQGKREKSDK DSKVSTSGS SVDSSKKTSE SKNVGSTGVA KIIISKHDGG SPSIKAKVTL QKPGESSGEG LRPQMASSKN YGSPLISGST PKHERGSPSH SKSPAYTPQN LDSESESGSS IAEKSYQNSP SSDDGIRPLP EYSTEKHKKH KKEKKVKDK DRDRDRDKDR DKKKSHSIKP

image



Background

Alternative Names	Mediator of RNA polymerase II transcription subunit 1 (Activator-recruited cofactor 205 kDa component) (ARC205) (Mediator complex subunit 1) (Peroxisome proliferator-activated receptor-binding protein) (PBP) (PPAR-binding protein) (Thyroid hormone receptor-associated protein complex 220 kDa component) (Trap220) (Thyroid receptor-interacting protein 2) (TR-
--------------------------	--

Product Name: Recombinant Human MED1 Protein
Catalog #: PEH4922



interacting protein 2) (TRIP-2) (Vitamin D receptor-interacting protein
complex component DRIP205) (p53 regulatory protein RB18A)

Note

For Research Use Only , Not for Diagnostic Use