


Product Name: Recombinant Human ZBTB6 Protein
Catalog #: PEH4925



Summary

Name	ZBTB6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15916
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFEQQGDVVL QKMNLLRQQN LFCDVSIYIN DTEFQGHKVI LAACSTFMRD QFLLTQSKHV RITILQSAEV GRKLLSCYT GALEVKRKEL LKYLTAASYL QMVHIVEKCT EALSKYLEID LSMKNNNQHT DLCQSSDPDV KNEDENSDKD CEIIEISEDS PVNIDFHVKE EESNALQSTV ESLTSERKEM KPELSTVDI GFKDNEICIL

image



Background

Alternative Names Zinc finger and BTB domain-containing protein 6 (Zinc finger protein 482)
(Zinc finger protein with interaction domain)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZBTB6 Protein
Catalog #: PEH4925

