


Product Name: Recombinant Human ZNF79 Protein
Catalog #: PEH4926



Summary

Name	ZNF79_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15937
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QALSEASFQD PCVEMPPGDS DHGTSBLEKS FNLRPVLSLQ QRVPVEARPR KCETHTEFSK NSEILKPHRA KPYACNECGK AFSYCSSLSQ HQKSHTGEKP YECSECGKAF SQSSSLIQHQ RIHTGEKPYK CSECGRAFSQ NANLTKHQRT

image



Background

Alternative Names Zinc finger protein 79 (ZNFpT7)

Note

For Research Use Only , Not for Diagnostic Use