


Product Name: Recombinant Human ZSCAN26 Protein
Catalog #: PEH4928



Summary

Name	ZSC26_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16670
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NSKGLGQEPL CKQFRQLRYE ETTGPREALS RLRELCQQWL QPETHTKEQI LELLVLEQFL IILPKELQAR VQEHHPEARE DVVVVLEDLQ LDLGETGQQD PDQPKKQKIL VEEMAPLKGV QEQQVRHECE VTKPEKEKGE ETRIENGKLI VVTDSCGRVE SSGKISEPME AHNEGSNLER HQAKPKEKIE YKCSEREQRF IQHLDLIEHA STHTGKKLCE

image



Background

Alternative Names Zinc finger and SCAN domain-containing protein 26 (Protein SRE-ZBP) (Zinc finger protein 187)

Note

Product Name: Recombinant Human ZSCAN26 Protein
Catalog #: PEH4928



For Research Use Only , Not for Diagnostic Use