


Product Name: Recombinant Human MAN2A1 Protein
Catalog #: PEH4929



Summary

Name	MA2A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16706
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLAENNEIIS NIRDVINLS ESVEDGPKSS QSNFSQGAGS HLLPSQLSLS VDTADCLFAS QSGSHNSDVQ MLDVYSLISF DNPDGGVWKQ GFDITYESNE

image



Background

Alternative Names	Alpha-mannosidase 2 (EC 3.2.1.114) (Golgi alpha-mannosidase II) (AMan II) (Man II) (Mannosidase alpha class 2A member 1) (Mannosyl-oligosaccharide 1,3-1,6-alpha-mannosidase)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use