


Product Name: Recombinant Human ZNF827 Protein
Catalog #: PEH4930



Summary

Name	ZN827_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q17R98
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HWYGNSETP SEASYGEVQE NYKLSLEDRI QEQSTSPDTS LGSTTPSSHT LELVALDSEV LRDSLQCQDH LSPGVSSLCD DDPGSNKPLS SNLRRLEAG SLKLDAAATA NGRVESPVNV GSNLSFSPPS HHAQQLSVLA RKLAEKQEQN

image



Background

Alternative Names Zinc finger protein 827

Note

For Research Use Only , Not for Diagnostic Use