


Product Name: Recombinant Human KRT24 Protein
Catalog #: PEH4936



Summary

Name	K1C24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2M2I5
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQTMQNLNDR LANYLDKVRA LEEANTDLEN KIKEWYDKYG PGSGDGGSGR DYSKYYSIIE DLRNQIIAAT VENAGIILHI DNARLAADDF RLKYENELCL RQSVEADING LRKVLDDLTM TRSDLEMQIE SFTEELAYLR KNHEEEMKMN QGSSGGEVTV EMNAAPGTDL TKLLNDMRAQ YEELAEQNRN EAEERFNKQS ASLQAQISTD AGAATSAKNE ITELKRTLQA LEIELQSQLA MKSSLEGTLA DTEAGYVAQL SEIQTQISAL EEEICQIWGE TKCQNAEYKQ LLDIKTRLEV EIETYRRLLD

image



Background

Alternative Names	Keratin, type I cytoskeletal 24 (Cytokeratin-24) (CK-24) (Keratin-24) (K24) (Type I keratin-24)
--------------------------	---

Product Name: Recombinant Human KRT24 Protein
Catalog #: PEH4936



Note

For Research Use Only , Not for Diagnostic Use