


Product Name: Recombinant Human ZNF181 Protein
Catalog #: PEH4938



Summary

Name	ZN181_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2M3W8
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVIIKVVVKQ SYEFSNSKKN LEYIEKLEGK HGSQVDHFRP AILTSRESPT ADSVYKYNIF RSTFHSKSTL SEPQKISAEG NSHKYDILKK NLPKKSIVKN EKVNGGKKLL NSNKSAAFS

image



SDS



WB

Background

Alternative Names Zinc finger protein 181 (HHZ181)

Note

For Research Use Only , Not for Diagnostic Use