


Product Name: Recombinant Human POLR1F Protein
Catalog #: PEH4944



Summary

Name	RPA43_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3B726
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IREQLDAELL RYSELLGVP IAYDNIKVVG ELGDIYDDQG HIHLNIEADF VIFCPEPGQK LMGIVNKVSS SHIGCLVHGC FNASIPKPEQ LSAEQWQTME INMGDELEFE VFRLSDAAG VFCIRGKLN I TSLQFKRSEV SEEV TENGTE EAAKKPKKKK KKKDPETYEV DSGTTKLADD ADDTPMEESA LQNTNNANGI WEEEPKKKKK KKKHQEVQDQ DPVFQGS DSS GYQSDHKKKK KKRKHSEEA E

image



Background

Alternative Names DNA-directed RNA polymerase I subunit RPA43 (DNA-directed RNA polymerase I subunit F) (Twist neighbor protein)

Note

Product Name: Recombinant Human POLR1F Protein
Catalog #: PEH4944



For Research Use Only , Not for Diagnostic Use