


Summary

Name	ZSC23_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3MJ62
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EMQEGLLAVK VKEEEEHSC GPESGLSRNN PHTREIFRRR FRQFCYQESP GPREALQRLQ ELCHQWLRPE MHTKEQILEL LVLEQFLTIL PEELQAWVRQ HRPVSGEEAV TVLEDLEREL DDPGEQVLSH AHEQEEFVKE KATPGAAQES SNDQFQTL EE QLGYNLREVC PVQEIDGKAG TWNVELAPKR EISQEVKSLI QVLGKQNGNI TQIPEYGDTC DREGRLKQR VSSSVERPYI CSECGKSFTQ NSILIEHQRT HTGEKPYECD

image



Background

Alternative Names	Zinc finger and SCAN domain-containing protein 23 (Zinc finger protein 390) (Zinc finger protein 453)
--------------------------	----------------------------------------------------------------------------------------------------------

Product Name: Recombinant Human ZSCAN23 Protein
Catalog #: PEH4946



Note

For Research Use Only , Not for Diagnostic Use