


Product Name: Recombinant Human GNPTAB Protein
Catalog #: PEH4948



Summary

Name	GNPTA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3T906
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NDEEFKMQIT VEVDTREGPK LNSTAQKGYE NLVSPITLLP EAEILFEDIP KEKRFPKFKR HDVNSTRRAQ EEVKIPLVNI SLLPKDAQLS LNTLDLQLEH GDITLKGYNL SKSALLRSFL MNSQHAKIKN QAIITDETND SLVAPQEKQV HKSILPNSLG VSERLQRLTF PAVSVKVNGH DQGQNPPLDL ETTARFRVET HTQKTIGGNV TKEKPPSLIV PLESQMTKEK KITGKEKENS RMEENAENHI

image



Background

Alternative Names	N-acetylglucosamine-1-phosphotransferase subunits alpha/beta (EC 2.7.8.17) (GlcNAc-1-phosphotransferase subunits alpha/beta) (Stealth protein GNPTAB) (UDP-N-acetylglucosamine-1-phosphotransferase subunits alpha/beta) [Cleaved into: N-acetylglucosamine-1-phosphotransferase
--------------------------	--

Product Name: Recombinant Human GNPTAB Protein
Catalog #: PEH4948



subunit alpha; N-acetylglucosamine-1-phosphotransferase subunit beta]

Note

For Research Use Only , Not for Diagnostic Use