


Product Name: Recombinant Human PARP15 Protein
Catalog #: PEH4950



Summary

Name	PAR15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q460N3
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTQGVVGTVS KPCFTAYEMK IGAITFQVAT GDIATEQVDV IVNSTARTFN RKSGVSRAIL EGAGQAVESE CAVLAAQPHR DFIITPGGCL KCKIIHVP GKDVRKTVTS VLEECEQRKY TSVSLPAIGT GNAGKNPITV ADNIIDAIVD FSSQHSTPSL KTVKVVIFQP ELLNIFYDSM KKRDLASLN

image



Background

Alternative Names	Protein mono-ADP-ribosyltransferase PARP15 (EC 2.4.2.-) (ADP-ribosyltransferase diphtheria toxin-like 7) (ARTD7) (B-aggressive lymphoma protein 3) (Poly [ADP-ribose] polymerase 15) (PARP-15)
--------------------------	--

Note

Product Name: Recombinant Human PARP15 Protein
Catalog #: PEH4950



For Research Use Only , Not for Diagnostic Use