


Product Name: Recombinant Human PLEKHN1 Protein
Catalog #: PEH4951



Summary

Name	PKHN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q494U1
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TRTSHLPES SVPSTVGCSS QHTPDQANS D RASIGRRRTE LRRSGSSRSP GSKARAEGRG PVTPLHLDLT QLHRLSLESS PDAPDHTSET SHSPLYADPY TPPATSHRRV TDVRGLEEFL SAMQSARGPT PSSPLPSVPV SVPASDPRSC SSGPAGPYLL SKKGALQSRA AQRHRGSAKD GGPQPPDAPQ

image



Background

Alternative Names	Pleckstrin homology domain-containing family N member 1 (PH domain-containing family N member 1) (Cardiolipin and phosphatidic acid-binding protein)
--------------------------	--

Note

Product Name: Recombinant Human PLEKHN1 Protein
Catalog #: PEH4951



For Research Use Only , Not for Diagnostic Use