


Product Name: Recombinant Human RFTN2 Protein
Catalog #: PEH4955



Summary

Name	RFTN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q52LD8
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLSGQSSESG IEEELHHESG QYQMEQNGSP TSSKSRKGEA SDNKLYTVFN AFDDDSTSWA YQEGILSMKV TRKGSVISTL DADWLELTF YKQGLSLID SFVFWETSKG EHLPKSLEGF FIYEEEGSGV PGSSRKGND A IVVEQWTVIE GCEIKTDYGP

image



Background

Alternative Names Raftlin-2 (Raft-linking protein 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RFTN2 Protein
Catalog #: PEH4955

