


Product Name: Recombinant Human TCEAL5 Protein
Catalog #: PEH4960



Summary

Name	TCAL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5H9L2
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGKREDEGEP GDEGQLEDEG NQEKQKSEGE EDKPQSEGKP ASQAKPESQP RAAEKRPAED YVPRKAKRKT DRGTDDSPKD SQEDLQERHL SSEEMMRECG DVSRAQEELR KKQKMGGFHW MQRDVQDPFA PRGQRGVRGV

image



Background

Alternative Names	Transcription elongation factor A protein-like 5 (TCEA-like protein 5) (Transcription elongation factor S-II protein-like 5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use