


Product Name: Recombinant Human HABP4 Protein
Catalog #: PEH4963



Summary

Name	HABP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JVS0
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAGAAMQESF GCVVANRFHQ LLDDESPFD ILREAERRRQ QQLQRKRRDE AAAAAGAGPR GGRSPAGASG HRAGAGGRRE SQKERKSLPA PVAQRPDSPG GGLQAPGQKR TPRRGEQQGW NDSRGPEGML

image



Background

Alternative Names Intracellular hyaluronan-binding protein 4 (IHABP-4) (IHABP4) (Hyaluronan-binding protein 4) (Ki-1/57 intracellular antigen)

Note

For Research Use Only , Not for Diagnostic Use