


Product Name: Recombinant Human ODR4 Protein
Catalog #: PEH4966



Summary

Name	ODR4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5SWX8
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IHSSKPKVKD AVQAVKRDIL NTVADRCEML FEDLLLNEIP EKKDSEKEFH VLPYRVFVPL PGSTVMLCDY KFDDESAEEI RDHFMEMLDH TIQIEDLEIA EETNTACMSS SMNSQASLDN TDDEQPKQPI KTTMLLKIQQ

image



Background

Alternative Names	Protein odr-4 homolog (hODR-4) (LAG1-interacting protein) (Transactivated by transforming growth factor beta protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use