


## Summary

<b>Name</b>	UT14C_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q5TAP6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	47 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	PLLTPMEKAS LQAMSLEEAK MHRAELQRRAR ALQSYEAKA RKEKKIKSKK YHKVVKKGKA KKALKEFEQL QKVNPTVALE EMEKIENARM MERMSLKHQN SGKWAKSKAI MAKYDLEARQ AMQEQLAKNK ELTQKLQVAS ESEEEGGTE VEELLVPHVA NEVQMNVDGP NPWMFRSCTS DTKEAATQED PEQVPELAAH EVSASEAER PVAEEILLR EFEERQSLRK RSELNQDAEP ASSQETKDSS SQEVLSELRA LSQKLKEKHQ SRKQKASSEG TVPQVQREEP APEEAPELLL QRSERVQTL ELEEKGKEDC FQNKELPRPV LEGQQSERTP NNRPDAPKEK KEKEQLINLQ NFLTTQSPSV RSLAVPTIIE ELEDEEERDQ

## image



## Background

**Alternative Names** U3 small nucleolar RNA-associated protein 14 homolog C

**Product Name: Recombinant Human UTP14C Protein**  
**Catalog #: PEH4970**

---



**Note**

For Research Use Only , Not for Diagnostic Use