



Product Name: Recombinant Human SOGA3 Protein
Catalog #: PEH4973



Summary

Name	SOGA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TF21
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHPQQQLQEQE EMQEEMEKLR EENETLKNEI DELRTEMDEM RDTFFFEEDAC QLQEMRHELE RANKNCRILQ YRLRKAERKR LRYAQTGEID GELLRSLEQD LKVAKDVSVR LHHELENVEE KRTTTEDENE KLRQQQLIEVE IAKQALQNEL EKMKELSLKR RGSKDLPKSE KKAQQTPTTE DNEDLKCQLQ FVKEEAALMR KKMAKIDKEK DRFEHELQKY RSFYGDLDSP LPKGEAGGPP STREAELKLR LRLVEEEANI LGRKIVELEV ENRGLKAELD DLRGDDFNGS ANPLMREQSE SLSELRQHLQ LVEDETELLR RNVADLEEQN KRITAE LNKY KYKSGGHDSA RHHDNAKTEA LQEELKAARL QINELSGKVM QLQYENRVLM SNMQRYDLAS HLGIRGSPRD SDAESDAGKK ESDDDSRPPH RKREGPIGGE

image



Background

Product Name: Recombinant Human SOGA3 Protein
Catalog #: PEH4973



Alternative Names

Protein SOGA3

Note

For Research Use Only , Not for Diagnostic Use