


Product Name: Recombinant Human ATF7IP2 Protein
Catalog #: PEH4977



Summary

Name	MCAF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5U623
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PESTVTSTVG DKKTDQMVFH LETNSNESH DKRQSDNILC SEDSGFVPVE KTPNLVNSVT SNNCADDILK TDECSRTSIS NCESADSTWQ SSLDTNNNSH YQKKRMFSEN EENVKRMKTS EQINENICVS LERQTAFLEQ VRHLIQQEY SINYELFDKK LKELNQRIGK TECRNKHEGI ADKLLAKIAK

image



Background

Alternative Names	Activating transcription factor 7-interacting protein 2 (ATF7-interacting protein 2) (MBD1-containing chromatin-associated factor 2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATF7IP2 Protein
Catalog #: PEH4977

