


Product Name: Recombinant Human EMC10 Protein
Catalog #: PEH4978



Summary

Name	EMC10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5UCC4
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAGTGARGAG AEGREGEACG TVGLLLEHSF EIDDSANFRK RGSLLWNQQD GTLSSLQRQL SEEERGLRD VAALNGLYRV RIPRRPGALD GLEAGGYVSS FVPACSLVES HLSDQLTLHV DVAGNVVGVS VVTHPGGCRG HEVEDVDLEL FNTSVQLQPP TTAPGPETAA FIERLEMEQA QKAKNPQEQK

image



Background

Alternative Names	ER membrane protein complex subunit 10 (Hematopoietic signal peptide-containing membrane domain-containing protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EMC10 Protein
Catalog #: PEH4978

