


Product Name: Recombinant Human KRT79 Protein
Catalog #: PEH4985



Summary

Name	K2C79_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5XKE5
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IDPEIQRVRT QEREQIKTLN NKFASFIDKV RFLEQQNKVL ETKWALLQEQ GQNLGVTRNN LEPLFEAYLG SMRSTLDRLQ SERGRLDSEL RNVQDLVEDF KNKYEDEINK HTAAENEFVV LKKDVDAAYM GRMDLHGKVG TLTQEIDFLQ QLYEMELSQV QTHVSNTNVV LSMDNRRNLD LDSIAEVKA QYELIAQRSR AEAEAWYQTK YEELQVTAGK HGDNLRDTKN EIAELTRTIQ RLQGEADAAK KQCQQLQTAI AEAEQRGELA LKDAQKKGD LDVALHQAKE DLTRLLRDYQ ELMNVKLALD VEIATYRKL ESEESRMSG

image



Background

Alternative Names Keratin, type II cytoskeletal 79 (Cytokeratin-79) (CK-79) (Keratin-6-like) (Keratin-6L) (Keratin-79) (K79) (Type-II keratin Kb38)

Product Name: Recombinant Human KRT79 Protein
Catalog #: PEH4985



Note

For Research Use Only , Not for Diagnostic Use