


Product Name: Recombinant Human ARL6IP4 Protein
Catalog #: PEH4988



Summary

Name	AR6P4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q66PJ3
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	序列

image



Background

Alternative Names	ADP-ribosylation factor-like protein 6-interacting protein 4 (ARL-6-interacting protein 4) (Aip-4) (HSP-975) (HSV1-binding protein) (SR-15) (SRp25) (SR-25) (Splicing factor SRp37)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use