


Product Name: Recombinant Human CREB3L3 Protein
Catalog #: PEH4989



Summary

Name	CR3L3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68CJ9
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ICA EKPADPV DLS PRCNLTV KDLLSGSSG DLQQHHLGAS YLLRPGAGHC QELVLTEDEK KLLAKEGITL PTQLPLTKYE ERVLKKIRRK IRNKQSAQES RKKKKKEYIDG LETRMSACTA QNQELQRKVL HLEKQNLSSL EQLKKLQAIV

image



Background

Alternative Names	Cyclic AMP-responsive element-binding protein 3-like protein 3 (cAMP-responsive element-binding protein 3-like protein 3) (Transcription factor CREB-H) [Cleaved into: Processed cyclic AMP-responsive element-binding protein 3-like protein 3]
--------------------------	--

Note

Product Name: Recombinant Human CREB3L3 Protein
Catalog #: PEH4989



For Research Use Only , Not for Diagnostic Use