


Summary

Name	ARID2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68CP9
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLSKNIPNHK TSNHVGNGEI SPMEPQGTLD ITQQDTAKGD QLERISNGPV LTLGGSSVSS IQEASNAATQ QFSGTDLLNG PLASSLNSDV PQQRPSVVVS PHSTTSVIQG HQIIAVPDSG SKVSHSPALS SDVRSTNGTA ECKTVKRP AETDRETVAGI PNKVGVRIVT ISDPNNAGCS ATMVAVPAGA DPSTVAKVAI ESAVQQKQQH PPTYVQNVVP QNTPMPPSPA

image



Background

Alternative Names	AT-rich interactive domain-containing protein 2 (ARID domain-containing protein 2) (BRG1-associated factor 200) (BAF200) (Zinc finger protein with activation potential) (Zipzap/p200)
--------------------------	--

Product Name: Recombinant Human ARID2 Protein
Catalog #: PEH4990



Note

For Research Use Only , Not for Diagnostic Use