


Product Name: Recombinant Human GTF2IRD2B Protein
Catalog #: PEH4996



Summary

Name	GTD2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6EKJ0
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAVVVQGLPE GVAFQHPENY DLATLKWILE NKAGISFIIN RPFLGPESQL GGPGMVTDAA RSIVSPSESC GPINVKTEPM EDGSLKAE AVSVKKESED PNYYQYNMQG SHPSSTSNEV IEMELPMEDS TPLVPSEEPN EDPEAEVKIE GNTNSSSVTN SAAGVEDLNI VQVTVPDNEK

image



Background

Alternative Names	General transcription factor II-I repeat domain-containing protein 2B (GTF2I repeat domain-containing protein 2B) (Transcription factor GTF2IRD2-beta)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GTF2IRD2B Protein
Catalog #: PEH4996

