


Product Name: Recombinant Human GRAMD4 Protein
Catalog #: PEH4997



Summary

Name	GRAM4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IC98
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGHKRDDFLD LAESPASDT ECSDEIPLKV PRTSPRDSEE LRDPAGPGTL IMATGVQDFN RTEFDRLNEI KGHLEIALLE KHFLQEELRK LREETNAEML RQELDRERQR RMELEQKVQE VLKARTEEQM AQQPPKGQAQ ASNGAERRSQ GLSSRLQKWF YERFGEYVED FRFQPEENTV ETEEPLSARR

image



Background

Alternative Names GRAM domain-containing protein 4 (Death-inducing protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GRAMD4 Protein
Catalog #: PEH4997

