

Product Name: Recombinant Human ZSCAN25 Protein
Catalog #: PEH5002



Summary

Name	ZSC25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NSZ9
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCHRQGEQEE TALCRGAWEP GIQLGPVEVK PEWGMPPGEG VQGPDPGTEE QLSQDPGDET RAFQEALPV LQAGPLPAV NPRDQEMAAG FFTAGSQGLG PFKDMALAFP EEEWRHVTPA QIDCFGEYVE PQDCRVSPGG GSKEKEAKPP QEDLKGALVA LTSERFGEAS LQGPGGLGRVC EQEPGGPAGS APGLPPPQHG AIPLPDEVKT HSSFWKPFQC PECGKGFSRS SNLVRHQRTTH EEKSYGCVCEC GKGFTLREYL MKHQRTHLGK RPYVCSECWK

image



Background

Alternative Names Zinc finger and SCAN domain-containing protein 25 (Zinc finger protein 498)

Product Name: Recombinant Human ZSCAN25 Protein
Catalog #: PEH5002



Note

For Research Use Only , Not for Diagnostic Use