


Summary

Name	SCMC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NUK1
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTAACQDAEQ PTRYETLFQA LDRNGDGVVD IGELQEGLRN LGIPLGQDAE EKIFTTGDVN KDGKLDFFEF MKYLKDHEKK MKLAFKSLDK NNDGKIEASE IVQSLQTLGL TISEQQAELI LQSIDVDGTM TVDWNEWRDY FLFNPVTDIE EIIRFWKHST GIDIGDSLTI PDEFTEDK

image



Background

Alternative Names	Calcium-binding mitochondrial carrier protein SCaMC-1 (Mitochondrial ATP-Mg/Pi carrier protein 1) (Mitochondrial Ca(2+)-dependent solute carrier protein 1) (Small calcium-binding mitochondrial carrier protein 1) (Solute carrier family 25 member 24)
--------------------------	--

Product Name: Recombinant Human SLC25A24 Protein
Catalog #: PEH5004



Note

For Research Use Only , Not for Diagnostic Use