


Summary

Name	PR40B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NWX9
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGSEDCDVLE ATQPLEQGFL QQLEEGPSSS GQHQPQQEEE ESKPEPERSG LSWSNREKAK QAFKELLRDK AVPSNASWEQ AMKMOVTDPR YSALPKLSEK KQAFNAYKAQ REKEEKKEAR LRAKEAKQTL QHFLEQHERM TSTTRYRRAE QTFGELEVWA VVPERDRKEV YDDVLFFLAK KEKEQAKQLR RRNIQALKSI LDGMSSVNFQ TTWSQAQQYL MDNPSFAQDH QLQNMDKEDA LICFEEHIRA LEREEEEERE

image



Background

Alternative Names	Pre-mRNA-processing factor 40 homolog B (Huntingtin yeast partner C) (Huntingtin-interacting protein C)
--------------------------	--

Product Name: Recombinant Human PRPF40B Protein
Catalog #: PEH5006



Note

For Research Use Only , Not for Diagnostic Use