

Product Name: Recombinant Human ZNF449 Protein
Catalog #: PEH5011



Summary

Name	ZN449_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6P9G9
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASLNQGSVFQ EYDTCDEVFR QRFRQFYRE AAGPHEAFNK LWELCCQWLK PKMRSKEQIL ELLVLEQFLT ILPTEIETWV REHCPENRER VVSLIEDLQR ELEIPEQQVD MHDMLLEELA PVGTAHIPPT MHLESPALQV MGPAQEAPVA EAWIPQAGPP ELNYGATGEC QNFLDPGYPL PKLDMNFSLE NREEPWVKEL QDSKEMKQLL DSKIGFEIGI ENEEDTSKQK KMETMYPFIV TLEGNALQGP ILQKDYVQLE NQWETPPEDL QTDLAKLVDQ QNPTLGETPE NSNLEEPLNP

image



Background

Alternative Names	Zinc finger protein 449 (Zinc finger and SCAN domain-containing protein 19)
--------------------------	---

Product Name: Recombinant Human ZNF449 Protein
Catalog #: PEH5011



Note

For Research Use Only , Not for Diagnostic Use