



Product Name: Recombinant Human TRAT1 Protein
Catalog #: PEH5012



Summary

Name	TRAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PIZ9
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SHYVEKQRQD KMSYSSDHT RVDEYYIEDT PIYGNLDDMI SEPMDENCYE QMKARPEKSV NKMQEATPSA QATNETQMCY ASLDHSVKGK RRPKRQKNTH FSDKDGDEQL HAIDASVSKT TLVDSFSPES QAVEENIHDD PIRLFGLIRA

image



Background

Alternative Names	T-cell receptor-associated transmembrane adapter 1 (T-cell receptor-interacting molecule) (TRIM) (pp29/30)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use