


Product Name: Recombinant Human FOXR2 Protein
Catalog #: PEH5013



Summary

Name	FOXR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PJQ5
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGVMKPPEMP QKRRSPDGD GPPCEPNLWM WVDPNILCPL GSQEAPKPSG KEDLTNISPF PQPPQKDEGS NCSEDKVVES LPSSSEQSP LQKQGIHSPS DFELTEEEAE EPDDNSLQSP

image



Background

Alternative Names Forkhead box protein R2 (Forkhead box protein N6)

Note

For Research Use Only , Not for Diagnostic Use