


Product Name: Recombinant Human ATAD2 Protein
Catalog #: PEH5015



Summary

Name	ATAD2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PL18
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTEETQDTSV DHNETGNTGE SSV EENEKQQ NASESKLELR NNSNTCNIEN ELEDSRKT TA CT ELRDKIAC NGDASSQII HISDENEGKE MCVLRMTRAR RSQVEQQQLI TVEKALAILS QPTPSLVVDH ERLKNLLKTV VKKSQNYNIF QLENLYAVIS QCIYRHRKDH DKTSLIQKME

image



Background

Alternative Names ATPase family AAA domain-containing protein 2 (EC 3.6.1.-) (AAA nuclear coregulator cancer-associated protein) (ANCCA)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ATAD2 Protein
Catalog #: PEH5015

