


Product Name: Recombinant Human GPR179 Protein
Catalog #: PEH5016



Summary

Name	GP179_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PRD1
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKALTRSRST YREKESVEES PEGQNSGTAG ESMGAPSRSP RLGRPKAVSK QAALIPSDDK ESLQNQQNAH TSRMLQVCQR EGSREQEDRG RRMTQGLGER KAERAGKTGL AMLRQVSRDK NIKQSKETPV GWQELPKAGL

image



Background

Alternative Names	Probable G-protein coupled receptor 179 (Probable G-protein coupled receptor 158-like 1) (GPR158-like)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use