


Product Name: Recombinant Human USP17L2 Protein
Catalog #: PEH5017



Summary

Name	U17L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6R6M4
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRATQGELKR DHPCLQAPEL DERLVERATQ ESTLDHWKFP QEQNKTKEPF NVRKVEGTLP PNVLVHQSK YKCGMKNHHP EQQSSLLNLS STTRTDQESV NTGTLASLQG RTRRSKGKKNK

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 17 (USP17) (EC 3.4.19.12) (Deubiquitinating enzyme 17-like protein 2) (Deubiquitinating protein 3) (DUB-3) (Ubiquitin carboxyl-terminal hydrolase 17-like protein 2) (Ubiquitin thioesterase 17-like protein 2) (Ubiquitin-specific-processing protease 17-like protein 2)
--------------------------	--

Product Name: Recombinant Human USP17L2 Protein
Catalog #: PEH5017



Note

For Research Use Only , Not for Diagnostic Use