


Product Name: Recombinant Human SAMD1 Protein
Catalog #: PEH5018



Summary

Name	SAMD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6SPF0
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLGALALPRG DRPGRAPPAA SARPSRSKRG GEERVLEKEE EEDDDEDEDE EDDVSEGSEV PESDRPAGA Q HHQLNGERGP QSAKERVKEW TPCGPHQGQD EGRGPAPGSG TRQVFSMAAM

image



Background

Alternative Names Atherin (Sterile alpha motif domain-containing protein 1) (SAM domain-containing protein 1)

Note

For Research Use Only , Not for Diagnostic Use