


Product Name: Recombinant Human ANKRD12 Protein
Catalog #: PEH5019



Summary

Name	ANR12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UB98
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFHLEFGEKS NAKIKDEDHS PTFENSCTL KKMDKEGKTL KKHKLKHKER EKEKHKEIE GEKEYKTKD SAKELQRSVE FDREFWKENF FKSDETEDLF LNMEHESLTL EKSKLEKNI KDDKSTKEKH VSKERNFKEE RDKIKKESEK SFREEKIKDL KEERENIPTD KDSEFTSLGM SAIEESIGLH LVEKEIDIEK QEKHIKESKE KPEKRSQIKE KDIEKMERT FEKEKKIKHE HKSEKDKLDL SECVDKIKEK DKLYSHHTEK CHKEGEKSKN TAAIKTDDR EKSREKMDRK HDKEKPEKER HLAESKEKHL MEKKNKQSDN SEYSKSEK GK NKEKDRELDK KEKSRDKESI

image



Background

Alternative Names	Ankyrin repeat domain-containing protein 12 (Ankyrin repeat-containing cofactor 2) (GAC-1 protein)
--------------------------	--

Product Name: Recombinant Human ANKRD12 Protein
Catalog #: PEH5019



Note

For Research Use Only , Not for Diagnostic Use