


Product Name: Recombinant Human ANKRD11 Protein
Catalog #: PEH5020



Summary

Name	ANR11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UB99
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLKEDKISKE K EKIFKEDKE K LKKEKVYRE DSAFDEYCNK NQFLENEDTK FSLSDQQRDR WFSDDLSDSSF DFKGEDSWDS PVTDYRDMKS DSVAKLILET VKEDSKERRR DSRAREKRDY REPFFRKKDR DYLDKNSEKR KEQTEKHKSV PGYLSEKDKK RRESAEAGRD RKDALESCKE RRDGRAKPPE AHREELKECG

image



Background

Alternative Names	Ankyrin repeat domain-containing protein 11 (Ankyrin repeat-containing cofactor 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ANKRD11 Protein
Catalog #: PEH5020

