


**Product Name: Recombinant Human FIP1L1 Protein**  
**Catalog #: PEH5021**



## Summary

<b>Name</b>	FIP1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6UN15
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	14 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LDAPGSINGV PLLEVDLDSF EDKPWRKPGA DLSDYFNYGF NEDTWKAYCE KQKRIRMGLE VIPVTSTTNK ITAEDCTMEV TPGAEIQDGR FNLFKVQQGR TGNSEKETAL

## image



## Background

**Alternative Names** Pre-mRNA 3'-end-processing factor FIP1 (hFip1) (FIP1-like 1 protein) (Factor interacting with PAP) (Rearranged in hypereosinophilia)

## Note

For Research Use Only , Not for Diagnostic Use