


Product Name: Recombinant Human ZNF66 Protein
Catalog #: PEH5029



Summary

Name	ZNF66_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZN08
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEKPYKCEEC GKAFSRSSIL TTHKIIHTGE KPYECEDCGK AFNRSSNLTK HKKIHTGEKP YKCEECGKAF KCSSILTTHK RIHTADKPYK CEECGKDFKY SSTLTRHKKI HTGGKPHKCN KCGKAFISSS

image



Background

Alternative Names Putative zinc finger protein 66 (Zinc finger protein 66 pseudogene)

Note

For Research Use Only , Not for Diagnostic Use