


Product Name: Recombinant Human BRNP3 Protein
Catalog #: PEH5038



Summary

Name	BRNP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q76B58
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CPSMDIQAME ENLLRITETW KAYNSDFEES DEFKLFMKRL PMNYFLNTST IMHLWTMDSN FQRRYEQLN SMKQLFLKAQ KIVHKLFSLs KRCHKQPLIS LPRQRTSTYW LTRIQSFLYC

image



Background

Alternative Names BMP/retinoic acid-inducible neural-specific protein 3 (DBCCR1-like protein 1)

Note

For Research Use Only , Not for Diagnostic Use