


Product Name: Recombinant Human CEP68 Protein
Catalog #: PEH5039



Summary

Name	CEP68_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q76N32
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAERSEPALS GLPPATMGSG DLLSGESQV EKTKLSSSEE FPQTLSPRT TTICSGHDAD TEDDPSLADL PQALDLSQQP HSSGLSCLSQ WKSVLSPGSA AQPSSCSISA SSTGSSLQGH QERAEPGRGS LAKVSSSLEP VVPQEPSSVV

image



Background

Alternative Names Centrosomal protein of 68 kDa (Cep68)

Note

For Research Use Only , Not for Diagnostic Use