


Product Name: Recombinant Human SUPT6H Protein
Catalog #: PEH5040



Summary

Name	SPT6H_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7KZ85
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IRVRPKKTTK KRVSRRSIFE MYEPSELESS HLTDQDNEIR ATDLPERFQL RSIPVKGAED DELEEEADWI YRNAFATPTI SLQESCDYLD RGQPASSFSR KGPSTIQKIK

image



Background

Alternative Names	Transcription elongation factor SPT6 (hSPT6) (Histone chaperone suppressor of Ty6) (Tat-cotransactivator 2 protein) (Tat-CT2 protein)
--------------------------	---------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use