


Summary

Name	GAK6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7LDI9
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKSQKETEGEGL HCEYVAEPVM AQSTQNVVDYN QLQEVYIPET LKLEGKGPPEL VGPSESKPRG TSPLPAGQVP VTLQPQKQVK ENKTQPPVAY QYWPPAELQY RPPPEQYGY PGMPPAPQGR APYPQPPTRR LNPTAPPSRQ GSKLHEIIDK SRKEGDTEAW QFPVTLEPMP PGEQAQEGEP PTVEARYKSF

image



Background

Alternative Names	Endogenous retrovirus group K member 6 Gag polyprotein (HERV-K(C7) Gag protein) (HERV-K(HML-2.HOM) Gag protein) (HERV-K108 Gag protein) (HERV-K_7p22.1 provirus ancestral Gag polyprotein) (Gag polyprotein)
--------------------------	--

Note

Product Name: Recombinant Human ERVK-6 Protein
Catalog #: PEH5043



For Research Use Only , Not for Diagnostic Use