


Product Name: Recombinant Human USP17L1 Protein
Catalog #: PEH5045



Summary

Name	U17L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7RTZ2
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRAKQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKTKEPF NVGKVEGTLP PNALVIHQSK YKCGMKNHHP EQQSSLLNLS STTRTDQESM NTGTLASLQG RTRRAKGKNK HSKRALLVCQ

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 17-like protein 1 (EC 3.4.19.12) (Deubiquitinating enzyme 17-like protein 1) (Ubiquitin thioesterase 17-like protein 1) (Ubiquitin-specific-processing protease 17-like protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human USP17L1 Protein
Catalog #: PEH5045

