


Product Name: Recombinant Human VPS35L Protein
Catalog #: PEH5047



Summary

Name	VP35L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z3J2
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRDRDDNSVV GSDFEPWTK RGEILARYTT TEKLSINLFM GSEKGKAGTA TLAMSEKVRT RLEELDDFEE GSQKELLNLT QQDYVNRIEE LNQSLKDAWA

image



Background

Alternative Names VPS35 endosomal protein-sorting factor-like (Esophageal cancer-associated protein)

Note

For Research Use Only , Not for Diagnostic Use